

# Dr. D. Y. Patil Pratishthan's DR. D. Y. PATIL INSTITUTE OF ENGINEERING, MANAGEMENT & RESEARCH

# Approved by A.I.C.T.E, New Delhi, Maharashtra State Government, Affiliated to Savitribai Phule Pune University

Sector No. 29, PCNTDA, Nigidi Pradhikaran, Akurdi, Pune 411044.

Phone: 020-27654470,Fax: 020-27656566

Website: www.dypiemr.ac.in Email: principal@dypiemr.ac.in

## **Experiment No:**

<u>Title</u>: - Measurement of the free space wavelength of the microwave (for TE 10 mode) with the help of the X-band microwave test bench and verify with its theoretical calculation.

#### **Objective:** -

To determine the frequency & wavelength in a rectangular waveguide working on TE10 mode.

## **Instruments & Components: -**

Gunn power supply

Gunn Oscillator

**Isolator** 

PIN modulator

Frequency meter

Slotted section

Tunable probe

Wave guide stand

SWR meter

Matched termination.

#### Theory: -

Mode represents in wave guides as either

TE m, n/TM m, n

Where,

TE – Transverse electric,

TM – Transverse magnetic

m – Number of half wavelength variation in x directionn- Number of half wavelength variation in y direction

$$\lambda g/2 = d1-d2$$

where d1 and d2 are the distance between two minima and maxima.

It is having highest cutoff frequency hence dominant mode. For dominant mode TE10 mode in rectangular wave guide  $\lambda \theta$ ,  $\lambda g$  and  $\lambda c$  are related below.

$$1/\lambda_0^2 = 1/\lambda_g^2 + 1/\lambda_c^2$$

Where

 $\lambda_0$  is free space wavelength

 $\lambda_s$  is guide wavelength

 $\lambda_c$  is cutoff wavelength

For TE10 mode,  $\lambda_c = 2a/m$ 

where m = 1 in TE10 mode and 'a' is broad dimension of waveguide.

The following Relationship can be proved

$$C = f\lambda$$

Where

 $C= 3 \times 10^8 \text{ m/s}$  is velocity of light

#### **Procedure: -**

1) Set up the components and equipments as shown in fig.

- 2) Set the variable attenuator at no attenuation position.
- 3) First connect the matched termination after slotted section.
- 4) Keep the control knob of Gunn power supply as shown.

Gunn bias knob : fully anti- clockwise direction PIN bias knob : fully anti- clockwise direction

PIN Mod frequency: mid position

Mode switch: Int. mode

Keep the control knob of SWR meter as shown.

Range dB: 50 dB Crystal: 200 ohm

Mode switch: Normal mode Gain (coarse & fine): mid position

SWR/dB: dB position

5) Set the micrometer of Gunn oscillator at 10 cm position.




Figure: Setup for study of frequency & wave length measurement.

- 6) Switch on the Gunn power supply, SWR meter and cooling fan.
- 7) Observe the Gunn diode current corresponding to the various voltages controlled by the Gunn bias knob through the LCD, don't exceed the bias voltage above 10.5 volts.
- 8) Turn the meter switch of power supply to beam voltage position and set beam voltage at 300V with help of beam voltage knob, current around 15 to 20mA.
- 9) Tune the probe for maximum deflection in SWR meter.
- 10)Tune the frequency meter to get a 'dip' minimum reading on SWR LCD display and note down the frequency directly from frequency meter. Now you can detune the DRF meter.
- 11) Move the tunable probe along with the slotted line to get the maximum reading in SWR meter.

- 12) Move the tunable probe to a minimum gain position record the probe position i.e. d1.
- 13) Move the probe to next minimum position and record the probe position again i.e. d2.
- 14) Calculate the guide wavelength as twice the distance between two successive minimum positions obtained as above.

$$\lambda g = 2(d1-d2)$$

15) Measure the wave-guide inner broad dimension 'a' which will be around 22.86mm for X band.

$$\lambda c = 2a$$

- Calculate the frequency by following equation:  $F = c/\lambda = c (1/\lambda_g^2 + 1/\lambda_c^2)^{1/2}$  Where  $c = 3 \times 108$  meter/sec. i.e. velocity of light.
- 17) Verify with frequency obtained by frequency meter.
- 18) Above experiment can be verified at different frequencies

# **Calculations:**

# **Conclusion:**